## **ERTELALSOP** THE FIRST NAME IN LIQUID FILTRATION<sup>™</sup>

## CONFIDENTIAL

## SPECIFICATION SHEET

<u>Product Code</u> :	Zeta-Pak® Lenticular Filter Cartridge	
<u>Grade Designation</u> :	12PFY7M504P06	
Materials:	Cellulose/Diatomaceous Earth	
<u>Wet Strength Resin</u> :	Non-Hazardous, listed as safe under 21CFR 176.170	
<u>Physical Properties</u> :		
	Surface:	Mottled
	Color:	White to Off-White
	Flowrate:	13-26 GPM per Pak® Average
	Retention:	1.0 micron nominal
	Configuration:	6 Cells

<u>Components:</u>

Gasket:

Teflon

7.0 sq. ft.

Support Members:

Filtration Area:

Polypropylene Core

- \* Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.
- <u>Note:</u> Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research

ERTELALSOP • 132 FLATBUSH AVENUE • KINGSTON • NY • 12401 PHONE 845.331.4552 • FAX 845.853.1526 • ERTELALSOP.COM