



# ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

*Confidential*

## **SPECIFICATION SHEET**

<b><u>Product Code:</u></b>	Zeta-Pak® Lenticular Filter Cartridge
<b><u>Grade Designation:</u></b>	16PFY6B6E910
<b><u>Materials:</u></b>	Cellulose / High Purity CelPure® Diatomaceous Earth
<b><u>Wet Strength Resin:</u></b>	Non-Hazardous, listed as safe under 21CFR 176.170
<b><u>Physical Properties:</u></b>	
<b>Surface:</b>	Mottled
<b>Color:</b>	White to Off-White
<b>Flowrate:</b>	135 LPM per Pak® Average
<b>Retention:</b>	0.80 micron   0.25 micron nominal
<b>Configuration:</b>	10 Cells
<b>Filtration Area:</b>	2.32 m <sup>2</sup>
<b><u>Components:</u></b>	
<b>Gasket:</b>	Silicone
<b>Support Members:</b>	Polypropylene

- \* Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.

**Note:** Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research