



# ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

*Confidential*

## **SPECIFICATION SHEET**

**Product Code:** Zeta-Pak® Lenticular Filter Cartridge

**Grade Designation:** 12PFY5M503LP16

**Materials:** Sulfate Treated Cellulose/CelPure® High Purity DE

**Wet Strength Resin:** Non-Hazardous, listed as safe under 21CFR 176.170

### **Physical Properties:**

***Surface:*** Mottled

***Color:*** White to Off-White

***Flowrate:*** 35-70 GPM per Pak® Average

***Retention:*** 1.0 micron nominal

***Configuration:*** 16 Cells

***Filtration Area:*** 18.5 sq. ft.

### **Components:**

***Gasket:*** Buna-N

***Support Members:*** Polypropylene

- \* Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.

**Note:** Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research

*Celpure is a registered trademark of Advanced Minerals, Inc.*