## **ERTELALSOP** THE FIRST NAME IN LIQUID FILTRATION<sup>™</sup>

## **SPECIFICATION SHEET**

<u>Product Code</u> :	Zeta-Pak® Lenticular Filter Cartridge	
<u>Grade Designation</u> :	12PFY5M5016	
<u>Materials</u> :	Cellulose/Perlite	
<u>Wet Strength Resin</u> :	Non-Hazardous, listed as safe under 21CFR 176.170	
<u>Physical Properties</u> :	Surface:	Mottled
	Color:	White to Off-White
	Flowrate:	15 to 20 GPM per Pak® Average
	Retention:	1.0 micron nominal
	Configuration:	16 Cells
	Filtration Area:	18.5 sq. ft.
<u>Components:</u>	Gasket:	Buna-N
	Support Members:	Polypropylene

- \* Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.
- *Note:* Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research