



ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

Confidential

SPECIFICATION SHEET

<u>Product Code:</u>	Zeta-Pak® Lenticular Filter Cartridge
<u>Grade Designation:</u>	12PFY5M703P16
<u>Materials:</u>	Cellulose/High Purity Celpure® DE
<u>Wet Strength Resin:</u>	Epichlorohydrin, Non-Hazardous, listed as safe under 21CFR 176.170
<u>Physical Properties:</u>	
Surface:	Mottled
Color:	White to Off-White
Flowrate:	15-20 GPM per Pak® Average
Retention:	0.50 micron nominal
Configuration:	16 Cells
Filtration Area:	18.5 sq. ft.
<u>Components:</u>	
Gasket:	Buna-N
Support Members:	Polypropylene

- * Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- * Retention based on retaining 98% of a given micron particle.
- * All components are listed in the CFR as generally recognized as safe for contact with food.

Note: Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research

Celpure is a registered trademark of Advanced Minerals, Inc.

ERTELALSOP • 132 FLATBUSH AVENUE • KINGSTON • NY • 12401
PHONE 845.331.4552 • FAX 845.853.1526 • ERTELALSOP.COM

SPECIFICATION SHEET