



# ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

*Confidential*

## **SPECIFICATION SHEET**

**Product Code:** Bio-Pak® Lenticular Filter Cartridge

**Grade Designation:** 12POX4M953P09

**Materials:** Cellulose/High Purity CelPure® DE

**Wet Strength Resin:** Non-Hazardous, listed as safe under 21CFR 176.170

**Physical Properties:**

<b>Surface:</b>	Mottled
<b>Color:</b>	White to Off-White
<b>Flowrate:</b>	7-15 GPM per Pak® Average
<b>Retention:</b>	0.25 micron nominal
<b>Configuration:</b>	9 Cells
<b>Filtration Area:</b>	10.4 sq. ft.

**Components:**

<b>Gasket:</b>	EPDM
<b>Support Members:</b>	Polypropylene

- \* Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.

**Note:** Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research