



ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

Confidential

SPECIFICATION SHEET

<u>Product Code:</u>	Zeta-Pak® Lenticular Filter Cartridge
<u>Grade Designation:</u>	12POX7M703P06
<u>Materials:</u>	Cellulose/High Purity CelPure® DE
<u>Wet Strength Resin:</u>	Epichlorohydrin, Non-Hazardous, listed as safe under 21CFR 176.170
<u>Physical Properties:</u>	
<i>Surface:</i>	Mottled
<i>Color:</i>	White to Off-White
<i>Flowrate:</i>	5-8 GPM per Pak® Average
<i>Retention:</i>	0.50 micron nominal
<i>Configuration:</i>	6 Cells
<i>Filtration Area:</i>	6.9 sq. ft.
<u>Components:</u>	
<i>Gasket/O-Ring:</i>	Teflon Encapsulated Silicone
<i>Support Members:</i>	Polypropylene

- * Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- * Retention based on retaining 98% of a given micron particle.
- * All components are listed in the CFR as generally recognized as safe for contact with food.

Note: Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research

Celpure® is a registered Trademark of Advanced Minerals