



ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

Confidential

SPECIFICATION SHEET

<u>Product Code:</u>	Zeta-Pak® Lenticular Filter Cartridge
<u>Grade Designation:</u>	12POY3M103P16
<u>Materials:</u>	Cellulose/High Purity Celpure® DE
<u>Wet Strength Resin:</u>	Epichlorohydrin, Non-Hazardous, listed as safe under 21CFR 176.170
<u>Physical Properties:</u>	
Surface:	Mottled
Color:	White to Off-White
Flowrate:	460 GPM per Pak® Average
Retention:	10 micron nominal
Configuration:	16 Cells
Filtration Area:	18.5 sq. ft.
<u>Components:</u>	
Gasket:	Viton
Support Members:	Polypropylene

- * Flowrate based on 70°F water flow through a 2 inch diameter filter at 10 PSI.
- * Retention based on retaining 98% of a given micron particle.
- * All components are listed in the CFR as generally recognized as safe for contact with food.
- * Celpure is a registered trademark of Advanced Minerals, Inc.

Note: Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research