



# ERTELALSOP

THE FIRST NAME IN LIQUID FILTRATION™

*Confidential*

## **SPECIFICATION SHEET**

**Product Code:** Zeta-Pak® Lenticular Filter Cartridge

**Grade Designation:** 16POY4M953P15NS

**Materials:** Cellulose/High Purity CelPure® DE

**Wet Strength Resin:** Epichlorohydrin, Non-Hazardous,  
listed as safe under 21CFR 176.170

### **Physical Properties:**

***Surface:*** Mottled

***Color:*** White to Off-White

***Configuration:*** 15 Cells

### **Performance Characteristics:**

***Flowrate:*** 30 GPM per Pak® Average  
114 LPM per Pak® Average

***Retention:*** 0.25 micron nominal

***Filtration Area:*** 37.5 sq. ft.  
3.5 sq. m.

### **Components:**

***Gasket:*** EPDM

***Support Members:*** Polypropylene

- \* Flowrate based on 70°F (21°C) water flow through a 2 inch (50mm) diameter filter at 10 PSI.
- \* Retention based on retaining 98% of a given micron particle.
- \* All components are listed in the CFR as generally recognized as safe for contact with food.
- \* CelPure® is a registered trademark of Advanced Minerals Corporation

**Note:** Manufactured in accordance with ErtelAlsop Drug Master File located at the United States Food and Drug Administration Center for Drug Evaluation and Research